

riva cem

레진 강화형 글래스아이오노머(RMGI) 루팅 시멘트

ionglass™



LIGHT TACK CURE (5초)로 빠르고 쉬운 잉여 시멘트 제거



PASTE 타입의 손쉬운 혼합 시스템






높은 본딩 강도



우수한 색조 안정성



지속적인 불소 방출



EVENT

riva cem 런칭 프로모션!

25.5g 8.5g (5.4ml) X 3 syringes (16.2ml)

정상가 ~~165,000원~~ → 할인가 **99,000원**

※ 프로모션 기간: 2018년 1월 ~ 3월

SDI

Your Smile. Our Vision.
www.sdi.com.au

구입 문의:

Key features




LIGHT TACK CURE로 빠르고 쉬운 잉여 시멘트 제거

고출력 LED 광중합기(460-480nm 출력)로 5초 광중합 후 가볍게 잉여시멘트를 제거할 수 있습니다.



우수한 색조 안정성

Colour changes in RMGICs stained with common drinks for 6 hours



PASTE 타입의 손쉬운 혼합 시스템



높은 본딩 강도

Shear bond strength to dentin (Mpa)



Shear bond strength to restorative (Mpa)




높은 압축 강도

Compressive strength (MPa)



지속적인 불소 방출

Cumulative fluoride release



Contents
 1 x riva cem 8.5g(5.4ml)
 1 x mixing pad
 1 x spatula

Technical Data	riva cem
Working time	1'45"
Setting time	4'30"
Shear bond strength (MPa) – Dentin	13.03 ± 3.40
Shear bond strength (MPa) – Enamel	9.99 ± 3.17
Compressive strength (MPa)	108

Indicated Material	type of restoration
PFM	Crowns & bridges
Resin	Crown, bridges, inlays, onlays & posts (fiber-reinforced)
Metals	Crown, bridges, inlays & onlays, orthodontic appliances
Ceramics (high strength)	Crowns & bridges
All-ceramic	Inlays

SDI

Your Smile. Our Vision.
 www.sdi.com.au

SDI CUBE 덴탈큐브

www.dentalcube.co.kr
 재료문의 : 02-2219-4567